

January 29, 2002

Amron International
759 West Fourth Avenue
Escondido
California – 92025

Product Name: Amron Oxygen Treatment Hood (Model –8891)

Test Name:

1. ISO Mem Elution - Passed
 2. ISO Kligman Maximization/44 animals - Passed
 3. ISO Intracutaneous/2 extracts - Passed
- The above tests were performed in compliance to and following the guidelines set forth in the current EN/ISO 10993. These tests meet the requirements of the Good Laboratory Practice and were conducted in compliance with U.S. Food and Drug Administration regulations set forth in 21 CFR, part 58.

Stacy Pritt DVM
Stacy Pritt, DVM
Study Director