

Internal and external measurements for the Avox Pressur-Vak II BIBS Mask.

Characteristic	01 – Small	02 – Medium	03 – Large	04 – XL	05-XS
A – inches (cm)	4.60 (11.7)	5.20 (13.2)	4.95 (12.6)	5.45 (13.8)	4.94 (12.5)
B – inches (cm)	4.91 (12.5)	4.94 (12.5)	5.13 (13)	5.05 (12.8)	4.65 (11.8)
C – inches (cm)	4.05 (10.3)	4.55 (11.6)	4.51 (11.5)	4.97 (12.6)	4.04 (12.6)
D – inches (cm)	4.20 (10.7)	4.82 (12.2)	4.70 (11.9)	4.62 (11.7)	4.20 (10.7)

A = Outside Width Measurement
 B = Outside Height Measurement
 C = Inside Height Measurement
 D = Inside Width Measurement

